

Prolystica® 2X Alkaline Cleaner




PERFORMANCE AND PROTECTION

Exceptional cleaning performance
and protection for automated
cleaning processes

NEW generation instrument cleaner

Prolystica® 2X Alkaline Cleaner

Soaking study to demonstrate the exceptional instrument protection provided by Prolystica 2X Alkaline Cleaner



With Prolystica



Without Prolystica

Features	Benefits
Sodium Hydroxide formulation with advanced surfactant system	Exceptional cleaning performance in the removal of blood, fat, mucus and other surgical soils from medical devices in all water types
Low available alkalinity	Broad range of material compatibility on instrument surfaces including anodized aluminium
Superior corrosion inhibitor system	Protects instruments from the harmful effects of water and water impurities
Enhanced Chelating and Sequestering system	Allows for efficient cleaning at high levels of water hardness
Free rinsing	Easily rinsed from instrument surfaces and doesn't require the use of an acid neutralizer
Concentrated Low Dose (1-4ml/L)	Less change over of product- more for less
Low foaming	Will not inhibit washer pump pressure and impede impingement
Environmentally Friendly	Has been formulated with biodegradable surfactants, chelating agents/sequestrants and corrosion inhibitors. Packaged in fully recyclable HDPE containers

STERIS environmental care

Prolystica Instrument Cleaning Chemistries have been formulated with biodegradable surfactants, chelating agents/sequestrants and corrosion inhibitors. Packaged in fully recyclable HDPE containers.



Please consult your STERIS representative for technical literature and product accessories.

Container sizes: 2 x 5 litres, 10 litres, and 20 litres.

STERIS OFFICES WORLDWIDE			
Benelux	32 2 523 2488	Japan	81 3 5521 2104
Canada	1 905 677 0863	Korea	82 2 517 1517
China	86 21 6137 1166	Latin America	1 800 884 9550
France	33 2 38 70 83 50	Malaysia	6 03 7954 9822
Germany	49 2203 890 6969	Singapore	65 68 41 7677
Greece	30 210 6800 848	Spain	34 91 658 5920
India	91 33 2367 5150	Switzerland	41 32 376 0200
Italy	39 022 13 034 1	United Kingdom	44 1256 840400

Document # M3341GB.2009-07, Rev. A
TPS Printed 09/2009, 1000

©2009 STERIS Corporation.
All rights reserved.

Technology to Prevent Infection and Contamination™

Printed on 100% recycled fibre.

STERIS®



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 ■ USA
440-354-2600 ■ 800-548-4873
www.steris.com

STERIS Limited
STERIS House
Jays Close
Viables
Basingstoke
Hampshire RG22 4AX ■ UK
+44 1256 840400 ■ 0800 252609
www.steris.com/international