

Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner





**PERFORMANCE, SPEED
AND PROTECTION**

Exceptional cleaning performance with speed
and protection for manual and automated
cleaning processes

NEW generation instrument cleaner

Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner



Features	Benefits
Dual Protease with advanced surfactant system	Exceptional cleaning performance with speed in the removal of blood, fat, mucus and other surgical soils from medical devices in all water types
Neutral pH	Broad range of material compatibility on instrument surfaces, including flexible endoscopes
Enhanced corrosion inhibitor system	Protects instruments from the harmful effects of water and water impurities
Formulated for manual & automated cleaning procedures	Single detergent solution
Concentrated Low Dose (1-4mL/L)	With a competitive cost per use
Low foaming	Will not inhibit washer pump pressure and impede impingement
Environmentally Friendly	Has been formulated with biodegradable surfactants, chelating agents/sequestrants and corrosion inhibitors. Packaged in fully recyclable HDPE containers
Free rinsing	Easily rinsed from instrument surfaces
Pleasant Fragrance	Deodorizes the smell of unpleasant surgical soils during the cleaning process

STERIS environmental care

Prolystica Instrument Cleaning Chemistries have been formulated with biodegradable surfactants, chelating agents/sequestrants and corrosion inhibitors. Packaged in fully recyclable HDPE containers.



Please consult your STERIS representative for technical literature and product accessories.

Container sizes: 6 x 1 litres, 2 x 5 litres, 10 litres and 20 litres



STERIS®



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 ■ USA
440-354-2600 ■ 800-548-4873
www.steris.com

STERIS Limited
STERIS House
Jays Close
Viables
Basingstoke
Hampshire RG22 4AX ■ UK
+44 1256 840400 ■ 0800 252609
www.steris.com/international

STERIS OFFICES WORLDWIDE	
Benelux	32 2 523 2488
Canada	1 905 677 0863
China	86 21 6137 1166
France	33 2 38 70 83 50
Germany	49 2203 890 6969
Greece	30 210 6800 848
India	91 33 2367 5150
Italy	39 022 13 034 1
Japan	81 3 5521 2104
Korea	82 2 517 1517
Latin America	1 800 884 9550
Malaysia	6 03 7954 9822
Singapore	65 68 41 7677
Spain	34 91 658 5920
Switzerland	41 32 376 0200
United Kingdom	44 1256 840400

Document # M3340GB.2009-07, Rev. A
TPS Printed 09/2009, 1000 ©2009 STERIS Corporation.
All rights reserved.

Technology to Prevent Infection and Contamination™

Printed on 100% recycled fibre.