

# RESERT™ XL HLD

## High Level Disinfectant



The safer, stronger, faster alternative to Glutaraldehyde and OPA from STERIS.

Ordering Information	Product Chamber	Quantity
Resert XL HLD	4450 AWEC	4 x 4 lt.
Verify Monitoring Strip for Resert	PCC030	2 x 60 each/bottle

# STERIS®



Powerful High Level Disinfection  
*including C.difficile spores*

Fast acting  
*8 minutes contact time at 20°C*

Long lasting Verify® Chemical  
Monitoring Strip for proof  
*21 days use*

Safe and easy to use  
*requires no special venting, measuring  
or mixing. No activation required*

Gentle on instruments  
*suitable for Endoscopes*

Aldehyde free  
*oxidative chemistry avoids protein fixing*

Non-toxic  
*formulated for single rinsing. Leaves no  
toxic residues or stains on devices*

 Safe  
*for staff, patients, instruments  
and the environment*



[www.steris.com/international](http://www.steris.com/international)

Tel: 0800 252609 Email: [UKHCsales@steris.com](mailto:UKHCsales@steris.com)